

VISION

We make it easier



GLOBAL A-SERIES™ Dental Microscope

All New Design Sleek, innovative style for the modern dental office

Intuitive Control The new AXIS™ Control System makes it easy to maneuver

Brighter Light With over 100,000 lux output, we offer the brightest LED light source available

Enhanced Comfort Updated ergonomic design to promote long term health benefits

GLOBAL

SURGICAL™ CORPORATION

We make it easier

Introducing the GLOBAL A-SERIES MICROSCOPE

Global is committed to delivering the best microscope experience available. The new A-Series™ microscope is designed for intuitive control while significantly improving visualization for earlier diagnosis, setting a new standard in dental microscopy.

- Available in 3, 4 and 6-steps of magnification
- AXIS™ Control System delivers magnification changer, handles and tension control in one place
- Updated binoculars with new locking diopter eyepieces
- Eyepieces available in 10x and 12.5x with softer eyecups
- Ergonomic handles with integral magnification changer
- Brightest LED light source available – 100,000 lux
- Three-stage filter assembly: clear, amber, green
- More mounting options than any other brand
- Lifetime warranty

MAGNIFICATION CHART*

Magnification (x)/Field of View (dia-mm)

A3	3.2 / 62	6.4 / 32	12.8 / 16			
A4	3.2 / 62	5.1 / 39	8.0 / 25	12.8 / 16		
A6	2.1 / 95	3.2 / 62	5.1 / 39	8.0 / 25	12.8 / 16	19.2 / 10

*Using 10x eyepieces, 250mm objective lens, 160mm binocular.
With selected components, this magnification range can be shifted up or down.



Floorstand



Wall Mount



High Wall Mount



Ceiling Mount