


CAMERA

ECO CAM X - DIGITAL VIDEO CAMERA 1-CCD

- Compact design
- High image quality
- Natural colour reproduction



CAMERA

ECO CAM X - DIGITAL VIDEO CAMERA 1-CCD

Technical Data

Camera head:	<ul style="list-style-type: none">• Aspect ratio: 4:3• Image sensor: 1/3" CCD• Immersible (using the attached cap)
System	<ul style="list-style-type: none">• Power supply: 100/240 VAC, 50-60 Hz• Fuses: 2 Fine fuses, 5x20 mm, delay acting 3,5 AH (I=1500A)• Power consumption: 15W• Image resolution: PAL 3x 752 x 582/ NTSC 3x 768 x 494• TV resolution: 470 lines• Sensivity: < 1 lux (1/1.4)• White balance automatically, with data hold• shutter speed PAL (NTSC): 1/50 (1/60) to 1/100.000 sec
Video out terminal	<ul style="list-style-type: none">• 1x Video (1Vpp/75 Ω), BNC• 1x S-Video (Y/C)• Mini-DIN-connector
Standard	<ul style="list-style-type: none">• Protection class: I• Application part: Type BF• IEC 601-1, EN 60601-1; classification acc. 93/42/EEC: I
Features	<ul style="list-style-type: none">• White balance• Window functions• Boost function• Anti-moiree-filter
Dimensions	<ul style="list-style-type: none">• Control unit: 176 x 72 x 202 mm (WxHxD); 1,6 kg
Delivery contains	<ul style="list-style-type: none">• CCU, camera head with 3 meter cable, mains cable• Y/C cable, BNC cable, user manual
Video system	<ul style="list-style-type: none">• PAL/NTSC

ecoCAM camera head



Weight: 80 g
Size: d= 29 mm; l= 55 mm

Without push buttons