

LIGHTSOURCE


HALOGEN HL 250

-
- Mechanical diaphragm
 - Softstart electronic
 - Exchangable lightguide adapter
 - Low heat conduction
-



LIGHTSOURCE

HALOGEN HL 250



Series	Description
HL 250	250 Watt, 2 lamps, 1 output
HL 250 DUO	250 Watt, 2 lamps, 2 outputs (no simultaneous use)

Technical Data	
Lamp	<ul style="list-style-type: none">• 250 W (1 or 2 Halogen reflector lamps), 24 V• 3600 K
Light guide connection	<ul style="list-style-type: none">• Storz, Wolf, Olympus, ACMI
System	<ul style="list-style-type: none">• Power supply: 110/115 VAC or 230/240 VAC switchable• Power consumption: 260 Watt
Standard	<ul style="list-style-type: none">• Protection class: I• Application part: Type BF• IEC 601-1, EN 60601-1; classification acc. 93/42/EEC: I
Dimensions	<ul style="list-style-type: none">• (WxHxD): 355 x 120 x 340 mm• Weight: 8 kg
Delivery contains	<ul style="list-style-type: none">• Light source, mains cable, user manual