

LIGHTSOURCE


LED 3000 - HIGH INTENSITY LIGHTSOURCE

-
- PHASER™ -Technology
 - Low noise, low heat conduction
 - Pure daylight quality
 - Waterproof touch control panel
-



LIGHTSOURCE

LED 3000 - HIGH INTENSITY LIGHTSOURCE



Series

Description

LED 3000-A	Storz adapter, automatic light intensity control
LED 3000-A	Wolf adapter, automatic light intensity control

Technical Data

Lamp	<ul style="list-style-type: none">• Daylight quality, 6000-6500 K• Lifetime: >10.000 hr
Light guide connection	<ul style="list-style-type: none">• Storz, Wolf adapter
System	<ul style="list-style-type: none">• Power supply: 100/240 VAC, 50-60 Hz• Fuses: 2 Fine fuses, 5x20 mm, delay 6,3 AT (I=1500A)• Power consumption: 300 W
Video In (opt.)	<ul style="list-style-type: none">• 1x Video, 1x S-Video (Y/C) (automatic Version only)
Standard	<ul style="list-style-type: none">• Protection class: I• Application part: Type BF• IEC 601-1, EN 60601-1; classification acc. 93/42/EEC: I
Dimensions	<ul style="list-style-type: none">• (WxHxD): 355 x 120 x 340 mm• Weight: 10,6 kg
Delivery contains	<ul style="list-style-type: none">• Light source, mains cable, user manual, BNC cable (optional)

Options

intensity displayed
as clear value

