

LIGHTSOURCE


XENON - BACKUP SERIES

- Low noise, low heat conduction
- Electronical light control via buttons
- Simple lamp replacement
- Waterproof touch control panel



LIGHTSOURCE

XENON - BACKUP SERIES




Series	Description
Xenon X 100-M-150	100 Watt, Storz adapter, manual light intensity control, backup lamp
Xenon X 180-M-150	180 Watt, Storz adapter, manual light intensity control, backup lamp
Xenon X 300-M-150	300 Watt, Storz adapter, manual light intensity control, backup lamp
Xenon X 100-MT-150	100 Watt, turret adapter, manual light intensity control, backup lamp
Xenon X 180-MT-150	180 Watt, turret adapter, manual light intensity control, backup lamp
Xenon X 300-MT-150	300 Watt, turret adapter, manual light intensity control, backup lamp
Xenon X 100-A-150	100 Watt, Storz adapter, automatic light intensity control, backup lamp
Xenon X 180-A-150	180 Watt, Storz adapter, automatic light intensity control, backup lamp
Xenon X 300-A-150	300 Watt, Storz adapter, automatic light intensity control, backup lamp
Xenon X 100-AT-150	100 Watt, turret adapter, automatic light intensity control, backup lamp
Xenon X 180-AT-150	180 Watt, turret adapter, automatic light intensity control, backup lamp
Xenon X 300-AT-150	300 Watt, turret adapter, automatic light intensity control, backup lamp

Technical Data

Lamp	<ul style="list-style-type: none">• Xenon arc lamp: 100W/ 180W/ 300W• Halogen backup lamp: 150W• Daylight quality, 5600 K• Lifetime: 500 hr (expected)
Light guide connection	<ul style="list-style-type: none">• Storz, Wolf, Olympus, or turret adapter
System	<ul style="list-style-type: none">• Power supply: 100-240 VAC, 50-60 Hz• Fuses: 2 Fine fuses, 5x20 mm, delay 6,3 AT (I=1500A)• Power consumption: 185W/ 280W/ 370W
Video In (opt.)	<ul style="list-style-type: none">• 1x Video, 1x S-Video (Y/C) (automatic Version only)
Standard	<ul style="list-style-type: none">• Protection class: I• Application part: Type BF• IEC 601-1, EN 60601-1; classification acc. 93/42/EEC: I
Dimensions	<ul style="list-style-type: none">• (WxHxD): 355 x 120 x 340 mm• Weight: 9,5kg/ 9,5kg/ 11kg
Delivery contains	<ul style="list-style-type: none">• Light source, mains cable, User manual, BNC cable (optional)

Options

automatic
light intensity control



turret connection
Storz, Wolf, Olympus, ACMI



lamp life/ intensity
displayed as clear value

