

# LIGHTSOURCE


## XENON - BASIC SERIES

- 
- Low noise, low heat conduction
  - Electronical light control via buttons
  - Simple lamp replacement
  - Waterproof touch control panel
- 



# LIGHTSOURCE

## XENON - BASIC SERIES




Series	Description
Xenon X 100-M	100 Watt, Storz adapter, manual light intensity control
Xenon X 180-M	180 Watt, Storz adapter, manual light intensity control
Xenon X 300-M	300 Watt, Storz adapter, manual light intensity control
Xenon X 100-MT	100 Watt, turret adapter, manual light intensity control
Xenon X 180-MT	180 Watt, turret adapter, manual light intensity control
Xenon X 300-MT	300 Watt, turret adapter, manual light intensity control
Xenon X 100-A	100 Watt, Storz adapter, automatic light intensity control
Xenon X 180-A	180 Watt, Storz adapter, automatic light intensity control
Xenon X 300-A	300 Watt, Storz adapter, automatic light intensity control
Xenon X 100-AT	100 Watt, turret adapter, automatic light intensity control
Xenon X 180-AT	180 Watt, turret adapter, automatic light intensity control
Xenon X 300-AT	300 Watt, turret adapter, automatic light intensity control

Technical Data	
Lamp	<ul style="list-style-type: none"> <li>· Xenon arc lamp: 100 W/ 180 W/ 300 W</li> <li>· Daylight quality, 5600 K</li> <li>· Lifetime: 500 hr (expected)</li> </ul>
Light guide connection	<ul style="list-style-type: none"> <li>· Storz, Wolf, Olympus, turret adapter</li> </ul>
System	<ul style="list-style-type: none"> <li>· Power supply: 100-240 V AC, 50-60 Hz</li> <li>· Fuses: 2 Fine fuses, 5x20 mm, delay 6,3 AT (I=1500A)</li> <li>· Power consumption: 185 W/ 280 W/ 370 W</li> </ul>
Video In (opt.)	<ul style="list-style-type: none"> <li>· 1x Video, 1x S-Video (Y/C) (automatic Version only)</li> </ul>
Standard	<ul style="list-style-type: none"> <li>· Protection class: I</li> <li>· Application part: Type BF</li> <li>· IEC 601-1, EN 60601-1; classification acc. 93/42/EEC: I</li> </ul>
Dimensions	<ul style="list-style-type: none"> <li>· (WxHxD): 355 x 120 x 340 mm</li> <li>· Weight: 9,0 kg/ 9,0 kg/ 10,5 kg</li> </ul>
Delivery contains	<ul style="list-style-type: none"> <li>· Light source, mains cable, user manual, BNC cable (optional)</li> </ul>

### Options

automatic light intensity control



turret connection Storz, Wolf, Olympus, ACMI



lamp life/ intensity displayed as clear value

