



ASKIR 30 SUCTION PUMP

NEW ASKIR 30 is a desk type electric suction unit suitable for nasal, oral and endotracheal aspiration of bodily fluids from adults and children. Indicated for tracheotomised patients, minor surgical applications and post operative therapy in home care. Main unit equipped with vacuum-meter (bar and kPa) and vacuum regulator. Complies with EN 60601-1-11 'for use in the home healthcare environment' and carries IP21 degree of protection.

Code	Description	Sell Price
W-3101	Askir 30 Suction Pump Complete	£320.00

Included Accessories:

Liquid Collection Jar with overflow valve system (different options)
 6 x 10 Silicone Tubes (autoclavable) + Conical Connector \varnothing 8-9-10mm
 Antibacterial & Hydrophobic Filter (single patient)
 CH20 Canula
 Power cord with Schuko plug

Specification:

Motor	Oilless and maintenance free piston pump
Power Feeding	230V-50Hz
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (Adjustable)	-0.88 Bar -80kPa -600 mmHG
Noise Level	60,5 dB
Max Free Air Flow Rate	40 l/min
Power Consumption	110 VA
Fuse	1 x F 1,6 A 250V
Duty Cycle	Non Stop Operation
IP Protection Code	IP21
Weight	3,6 KG
Size	35 x 21 x 18 cm
Warranty	2 Years
Place of Manufacture	Italy

All Prices are Subject to VAT and Delivery

DP Medical Systems Ltd
 15 a Oakcroft Road, Chessington, Surrey KT9 1RH
 Tel: 020 8391 4455 email: sales@dpmedicals.com www.dpmedicals.com