



Examination chairs for gynaecology



- Ergonomic and modern design
- Polished high quality stainless steel edition
- Optimal to clean and disinfect
- User friendly foot switch conception
- Additional hand switch with coloured display (400550/570)
- Accessible for handicapped and elderly people
- 12 memory positions / 4 user accounts
- Time saving by simultaneous drive of engines (memory-positions)
- Various colouring possible
- Small footprint by compact construction
- Maximum permissible overall load of 200 kg







Examination and treatment chair

Art.-no. 400520

With the following electrical functions:

- Height adjustment
- Adjustment of pelvis section
- Adjustment of back section

Shown accessories:

Art.-no. 40023 Neck upholstery

Art.-no. 81726_1 Footrests "Comfort" incl. artificial leather cover



Examination and treatment chair Art.-no. 400550

With the following electrical functions:

- Height adjustment
- Adjustment of pelvis section
- Adjustment of back section
- Hand switch
- Memory function with three programmable positions for four operators

Shown accessories:

Art.-no. 40023 Neck upholstery

Art.-no. 40000 Powder coated base according to RAL colour chart

Art.-no. 42000 Column cover made of stainless steel

Art.-no. 81717 Leg holders according to Goepel incl. artificial leather covering





Examination and treatment chair

Art.-no. 400570

With the following electrical functions:

- Height adjustment
- Adjustment of pelvis section
- Adjustment of back section
- Adjustment of leg rests
- Hand switch
- Memory function with three programmable positions for four operators

Shown accessories:

Art.-no. 40023 Neck upholstery

Art.-no. 81726 2 Footrests incl. artificial leather cover



Basic Features

- Aluminium base and column cover powder coated, white (RAL 9010)
- Seat part and back rest frame made of polished stainless steel
- Electric movements controlled by easy operated cable foot switch
- Electric movements controlled by easy operated hand switch with colour display (Model 400550, 400570)
- Floor levelling screws for secure standing
- Upholstery with ergonomic form, 50 mm thickness
- Head and back upholstery prepared for neck upholstery with magnetic fixation
- Plastic covered backsides of upholstery parts
- Paper roll holder for paper rolls (up to 500 mm) at back rest
- Paper rollable over back rest and seat part, or just over seat part
- Mounting device for side rails (25 x 10 mm), hand grips as well as a fixation device for leg holder system acc. to Goepel
- Extendable, inclinable, detachable white plastic basin
- Safe working load: 200 kg
- Weight of the chair: approx. 110 kg
- Floor space required: 860 x 550 mm
- Power supply: 110/230 V / 50-60Hz , incl. 3 m cable







Models

Comfort and Elegance

			200		
		Range	400520	400550	400570
\$	Height adjustment	550 - 900 mm	•	•	•
	Height adjustment by means of raising pelvis part	230 mm	5	50 - 1130 mm	
	Adjustment of pelvis section	+20°	•	•	•
*	Adjustment of back section	+50°	٠	٠	•
	Leg holder adjustment	30°			•
Μ	Memory (12 positions)			•	•
	Hand switch		•	•	•
ľ	Maximum permissible overall load		200 kg	200 kg	200 kg
			Electrical	Manual	 Optional

Upholstery

- upholstery parts

- Seamlessly manufactured
- and to disinfect
- any extra charge
- neck rest upholstery







Functionality and Design





Functionality and Design

Powered adjustment of height from 550 up to 900 mm

Powered adjustment of pelvis section up to + 20°

Powered adjustment of back section up to +50°



Powered adjustment of leg holder system up to 30°



Control concept

Learn more about the advantages of the medifa 4000 control concept using the model 400570 as an example.

- Only a few keys ensure simple and intuitive operation
- Control concept can be learned quickly thanks to the same symbols and colour scheme of the control units
- Colour display shows motion performed and selected user level (400550 and 400570 models only)
- Time saving due to synchronous running of the motors in memory mode
- Four user levels, each with three memory positions - ideal for multiple users
- All movements can be immediately stopped and executed again
- All three memory positions can also be activated directly via foot switch

Illuminated colour display shows the current position, of the 4 available user modes, of the 3 memory positions 	
Double function key for • preselection of the several functions • preselection of memory positions M1 to M3	
LED display shows selected function	Ruction
Positioning keys for all functions	
User selection key (4 users)	
Activating a saved memory position	









Ascent and descent position

14



Gynaecological treatment position

Ultrasonic position

Technical Information







1.	Head upholstery	6.	Power cable outlet
2.	Back upholstery	7.	Socket at the base (accessory 40677)
3.	Pelvis upholstery	8.	Leg plate (accessory 400076)
4.	Leg plate (accessory 400066)	9.	Colposcope holder, right or left (Accessory 82347_System)
5.	Potential equalisation (accessory 40678)	10.	Cable outlet for floor connection





R

Leg Holder Systems



Art.-no. 81726_1

Footrests Comfort (pair), made of stainless steel, with integrated hand grip function and artificial leather cover, for attachment at the chair's toothing elements

- Stepless, gas spring supported onehand height adjustment
- Pre-adjustment of the footrests by toothing elements at the chair
- Swivelling
- Grey foamed PUR footrests
- Not in combination with items 81072 and 82347_System







Art.-no. 81726_2

Footrests (pair), made of stainless steel, with integrated hand grip function and artificial leather cover, for attachment at the chair's toothing elements

- Adjustment of footrests by toothing elements at the chair
- Swivelling
- Grey foamed PUR footrests
- Suitable for electric adjustment (400570)

Leg Holder Systems



Art.-no. 81717

Leg holders according to Goepel (pair), made of stainless steel, with integrated hand grip function and artificial leather cover, for attachment at the chair's toothing elements (only for types 400520, 400550, 400570)

- Adjustment of leg holders by toothing elements at the chair
- Swivelling
- Grey foamed PUR footrests
- Stepless adjustment of the knee bolsters by ball joint
- Suitable for electric adjustment (400570)





Art.-no. 81705

Leg holder according to Goepel Comfort with grey foamed PUR knee bolsters (pair)

Art.-no. 85223

Side rails at pelvis section for accessory attachment (pair), I = 320 mm, only required for items 81705, 81700 and 81650

Art.-no. 81500

Clamps made of stainless steel, rotatable, for side rail attachment (pair)

Leg Holder Systems



Art.-no. 81700 Leg holders according to Goepel Classic with black foamed PUR knee bolsters (pair)

Art.-no. 81710 Additional charge for artificial leather cover of knee bolsters (only for item 81700) in upholstery colour





Art-no. 81650

Heel supports (pair) for fixation at side rail, made of stainless steel, without clamps





Art.-no. 400076

Leg plate for short horizontal positioning of the patient

- Complete insertion of the leg plate below the seat part
- Enlargement of the chair's lying surface up to 1670 mm
- Not in combination with items 400066, 91706 and 91767

Art.-no. 400066

Leg plate for horizontal positioning of the patient

- Leg plate adaption by insertion below the seat part
- Enlargement of the chair's lying surface up to 1950 mm
- Not in combination with item 400076

Gynaecological Accessories

Standard Accessories





Art.-no. 82347_System

Colposcope holder for systems Olympus, Kaps, Zeiss, Leisegang and Atmos, adaptable to the left or to the right, not in combination with item 81726 1

Art.-no. 81072

Gyn treatment set with extendible arm, swivelling, adaptable to the left or to the right, not in combination with items 81726 1 and 91716, incl. additional accessories: 1 stainless steel storing-peel for instruments, 3 glass cups with lids, 3 glass bottles with lids, 4 instrument holster, 1 object slide

Art.-no. 81051

Stainless steel board with extendible arm, swivelling, adaptable to the left or to the right, dimensions: L x W = 413 x 398 mm

Art.-no. 81047 200

Instrument drop made of stainless steel for fixation at column of chair, extendable by means of telescopic rail (Tip-on), with stainless steel bowl (12.5 litres) and stainless steel basket for specula, not in combination with item 40018

Dimensions (Stainless steel bowl): L x W x H = 300 x 235 x 200 mm Dimensions (Stainless steel basket): L x W x H = 230 x 175 x 140 mm







Art.-no. 40677 Socket at the base, 2 pieces (cannot be retrofitted)

Art.-no. 40678 Potential equalisation socket at base (cannot be retrofitted)

Art.-no. 40018 Doctor's foot rest made of stainless steel, not in combination with items 81047 200, 91762 and 91763

Art.-no. 42000 Column cover made of polished stainless steel

Art.-no. 40000 Additional charge for powder coated base according to RAL colour chart

Art.-no. 41000 Additional charge for powder coated cover of the column according to RAL colour chart

Standard Accessories





Mobility for examination chair, consisting of 4 double leading castors, ø =75mm, lockable, electrically conductive, not for patient transfer

Art.-no. 40046

Hand switch for control of height adjustment, pelvis section adjustment and back section adjustment (only for item 400520, cannot be retrofitted)

Art.-no. 40020 PVC-protective cover for pelvis part upholstery

Art.-no. 85222 Side rails at back section for accessory attachment (pair), I = 600 mm

Art.-no. 85223 Side rails at pelvis section for accessory attachment (pair), I = 320 mm, only required for items 81705, 81700 and 81650









Art.-no. 81043 Storing bowl made of stainless steel, dimensions: L x W x H: 325 x 265 x 65 mm

Art.-no. 81044 same as 81043, but with dimensions (L x W x H): 325 x 176 x 65 mm

Art.-no. 28530 Hand grips, white RAL 9010 (pair)

Art.-no. 28531 Additional charge for artificial leather cover of handgrips in colour of upholstery (pair)



Art.-no. 61211

Arm rest, rotatable and tiltable by means of ball joint, height adjustable, horizontally slidable, upholstery electrically conductive, incl. clamp and fixture strap, dimensions (upholstery): $L \times W \times H = 400 \times 115 \times 40$ mm

Art.-no. 40220

Arm rest, swivelling, incl. artificial leather covered upholstery in colour of choice, dimensions (upholstery): L x W = 600 x 165 mm

Standard Accessories



Art.-no. 40023 Neck upholstery, incl. magnetic fixation

Art.-no. 81805 Infusion rod for side rail attachment, stainless steel, incl. rotatable clamp





Art.-no. 40070

Wall holder made of stainless steel, for attachment of footrests, leg holders, leg plates and further accessories to fix at side rail (25 x 10 mm), dimensions (L x W x H): 300 x 235 x 200 mm



Art.-no. 82349 1 Holder for examination luminaire, adaptable on the right side

Art.-no. 82349 2 Holder for examination luminaire, adaptable on the left side

Art.-no. 82349_3 Examination luminaire holder for examination chair with colposcope holder, adaptable at the right or at the left side



Art.-no. 82348

Examination luminaire for attachment at holder 82349 1, 82349 2 or 82349 3

- LED technology (50'000 lx / 0.5 m)
- Illuminated field 180 mm
- Precise colour rendering Ra >93
- Light colours 3300 K, 3800 K or 4400 K
- Illumination strength dimmable in five stages
- Operating unit on light head
- Stepless positioning through flexible arm
- Slight thermal radiation in light

Additional Furniture

medifa

The medifa healthcare group unites competence concerning operating rooms under one roof.

Art.-no. 206020

Mobile stool with polished five-arm aluminium base, double swivel castors ø 50 mm, height adjustment with gas spring support from 420 to 550 mm via release lever under seat, upholstery made of artificial leather covered foam according to colour chart, dimensions: seat: ø 360 mm, upholstery: 60 mm, base: ø 600 mm

Art.-no. 206045

Mobile stool with saddle seat and base ring, polished five-arm aluminium base, double swivel castors ø 50 mm, height adjustment with gas spring support from 420 to 550 mm via release lever under seat, upholstery made of artificial leather covered foam according to colour chart, dimensions: seat: ø 400 mm. upholstery: 60 mm, base: ø 600 mm

medifa is an owner-run and internationally active complete solution provider in the medical technology.

The product and service portfolio covers all operating room requirements from modular wall, door, ceiling and cabinet systems (RooSy® Room Systems) for the turnkey operating room, as well as mobile operating tables, operating table accessories and operating lights, right up to examination chairs and medical

Products

- Modular operating room systems: walls, doors, ceilings and cabinets (RooSy)
- Mobile operating tables and accessories
- Operating lights
- Examination and treatment chairs
- Medical furniture and trolleys
- Customized medical solutions

furniture. In addition, one of the largest German OEM suppliers in the medical technology sector also belongs to medifa.

medifa GmbH & Co. KG (Finnentrop), medifa hygienic rooms GmbH (Ötigheim) and medifa metall und medizintechnik GmbH in Rastatt are all part of this group of companies. The medifa healthcare group employs approximately 350 employees at three locations within Germany. The planning, engineering and production

medifa worldwide





take place exclusively in Germany, and with over 90% added value in the group of companies, the group guarantees fast reaction times.

The slogan "we care" means for us that we consistently orientate ourselves according to the requirements of our customers. Customers all over the world can rely on top-class consultation, processing and quality from a single source, made in Germany according to certified national and international quality standards.



Made in Germany







Version 07-20-EN Subject to modifications due to technical development. This document is intended to provide information to an international audience outside of the US.

www.medifa.com



medifa GmbH & Co. KG

Industriestraße 5 57413 Finnentrop Germany fon +49 2721 71 77 -0 fax +49 2721 71 77 -255 info@medifa.com

