

DP
MEDICAL

VISION IN HEALTHCARE

Specialist equipment suppliers for FME and SARC

MediScan HD Patient Imaging System


- MediScan software provides a safe and secure pathway for recording, storing and encrypting examinations for forensic analysis and archive purposes.
- Our MediScan software can be supplied on a number of different hardware; from laptops to high definition touchscreen panel PC's for quick and easy use in OPD, Theatre and examination rooms.
- As well as videos and images, the software facilitates the capture of information that can be burned to encrypted DVD. The datasets are customisable for individual use, and can be used for admin purposes.

MediScan[®]

KAPS Vico S HD Video Colposcope

- Offering Full HD quality images and video, the KAPS Vico S HD Video Colposcope is ideal for Forensic Medical Examinations.
- The Colposcope is mounted on a mobile stand; and the spring arm allows for even easier positioning of the optics for unobstructed views of the patient.
- Featuring an array of LED bulbs, the illumination enables brilliant images with enhanced image detail.





MediScan Compact Video Colposcope



- Handheld Video Colposcope designed to be used with our MediScan software on a laptop.
- This versatile device can be used mounted on a mobile stand or alternatively as a hand held device by the clinician.
- Offering Full HD image and video quality; with high performance LED illumination.

Duolys Patient Couch

- An ideal couch for a multi-purpose clinical space / room.
- Motorised functions for height, and a central section that can be raised to transform the couch into a gynaecological & forensic consultation chair.
- Integrated retractable stirrups and smooth electrical movement provide the user and the patient with great comfort.

