



**enterMed**

**New  
Interview III**

Cordless examination headlight



**Clear, intense  
LED light**



**Universal charger  
'Micro-USB'**



**Comfortable  
& hygienic**

- Micro-USB-b
- Cordless, rechargeable design
- Lightweight (325 g)
- Focusable light beam
- Multipurpose
- Flexible lamp housing
- High battery capacity

## Enterview III technical specifications

Battery capacity:	Approx. 3 hours when fully charged
Weight:	325 grams
Spot diameter (distance 25 cm):	25 - 100mm
Light hue:	5650K
Charger connector:	Micro USB-B with status LED
Increase operation time:	Power bank 5V, min. 1A ("not included")
Luminous flux:	235 lumen
Collimating system:	Fresnel lens
Battery type:	Li-Ion cell, 2200 mAh with NTC protection
USB battery charger included.	

### Product codes:

**670600 Enterview III cordless headlight LED-type,  
including USB adaptor**

Entermed B.V.  
Edisonstraat 22  
8912 AW Leeuwarden  
The Netherlands  
Tel: +31 58 298 8888  
Email: [info@entermed.com](mailto:info@entermed.com)  
Internet: [www.entermed.com](http://www.entermed.com)

We continuously strive to improve and update our products and thus reserve the right to make alterations without prior notice.