

TDCL LED Medical Headlight



Usage

ENT, dental surgery, plastic surgery, orthopaedics, neurosurgery, cardiology, ophthalmology, gynaecology, veterinary.

Specifications

Voltage	AC90-240V
Lamp Power	LED 3W
Color temperature	5500±500K
Lamp life	50000h
Illumination	40000lx
Charging time	5h
Power supply time	4h/2h
N.W.	280g+80g
Package size	205*190*127mm

£595 plus VAT

Prices subject to delivery and VAT